

#### FEATURES AND BENEFITS

6.5L ● Dosing Capacity

#### Retaining Bolts on lid

- Avoids misplacing bolts
- Allows for dosing directly into the unit as well as via the tundish

## 2 options for both inlet & outlet connections

- Easy installation
- Saves time and costs





#### **Isolation Valves**

- Under AAV prevents leaking
- All stainless steel to enable use on chilled systems

### Single Large Magnetic Rod

- Large surface area for maximum ferrous capture
- Easy Cleaning Nickel plated avoids rusting and demagnetisation



#### Mesh Basket Filter

- 100 litres per minute flow rate
- Washable and reusable (No consumables i.e cartridge)
- Saves time and cost





# SS Magnetic Filter



Disclaimer: The information within this document is believed to be correct at the time of publication; however, the document is for guideline use only. For complete accuracy, always check the product with a CORE representative. Missing information was either not available or disclosed. It is your responsibility that any product meets the necessary requirements. Any reliance placed upon this information will be totally at the user's risk.







www.sbscore.co.uk

Order online

today



V4.